
# USB PEDOMETER

Quickstart Guide  
Guide d'utilisation rapide  
Schnellstart-anleitung

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### CAUTION:

- Before starting any exercise program, we recommend consulting your physician for a physical exam and to discuss your program.
- The pedometer is not water-resistant. Do not submerge it in liquid or expose it to wet conditions.
- Avoid rough use or severe impact to the pedometer.

### FEATURES

1. **MODE/DOWN Button** – Toggles between the available modes: (A) Time Mode, (B) Step Mode, and (C) Steps per Minute Mode.
2. **UP Button** – Increases the value of the current setting or toggles through available options.
3. **SET Button** – Saves the current setting.
4. **REC Button** – Toggles through the readings (steps, calories, distance, etc.) stored in the pedometer.
5. **LCD** – Displays information related to the pedometer's operation.

### ACTIVATE THE PEDOMETER:

Before using the pedometer for the first time, remove the clear tape from the battery compartment on the back. The pedometer will turn on (and remain on for the life of the battery, which can last four months to one year).

### SET YOUR PERSONAL DATA

#### UNITS OF MEASUREMENT

1. Press MODE to toggle through the available modes until Step Mode is shown.
  2. Press and hold SET until the LCD starts to flash (~10 seconds).
  3. Use the UP or MODE/DOWN buttons to toggle between "KM" (kilometers) and "ML" (miles).
  4. Press SET to save the setting.
- Note:** If you select "KM," your height will be measured in "CM" (centimeters), and your weight will be measured in "KG" (kilograms). If you select "ML," your height will be measured in "IN" (inches), and your weight will be measured in "LB" (pounds).

#### TIME, ALARM, AND USER INFORMATION

1. Press MODE to toggle through the available modes until Time Mode is shown.
2. Press and hold SET until the Hour setting starts to flash (~3 seconds).
3. Use the UP or MODE/DOWN buttons to increase/decrease the hour. Press SET to save the setting.
4. Repeat Step #3 for the Minute setting.
5. The Alarm setting will be shown with its hour number flashing. Use the UP or MODE/DOWN buttons to increase/decrease the hour. Press SET to save the setting.
6. Repeat Steps #5 for the Minute setting.
7. To turn the Alarm on/off, use the UP or MODE/DOWN buttons. (The "bell" icon indicates the Alarm is on.) Press SET to save the setting.
8. Use the UP or MODE/DOWN buttons to toggle between "12" (a 12-hour time system) or "24" (a 24-hour time system). Press SET to save the setting.
9. The Year number will be flashing. Use the UP or MODE/DOWN buttons to increase/decrease the Year. Press SET to save the setting.
10. Repeat Step #9 for the Month, then the Date.
11. Your Height setting will flash. Use the UP or MODE/DOWN buttons to increase/decrease your Height. Press SET to save the setting. (**Note:** If you are planning on running, see the *RUNNING* section below to set your height properly.)
12. Repeat Step #11 for your Weight, then your Age.
13. An icon representing your Gender will flash. Use the UP or MODE/DOWN buttons to toggle between the "male" and "female" icons. Press SET to save the setting.
14. The largest two decimal places of your target calorie goal will flash. Use the UP or MODE/DOWN buttons to increase/decrease it. Press SET to save the setting. (**Note:** "KCAL" stands for "kilocalories," which is actually the same as "food calories" (the calories you monitor in your diet).)
15. Repeat Step #14 for the smallest two decimal places of your calorie goal. Press SET to complete the setup.

\* If the Alarm is on, it will sound whether or not you have reached your calorie goal. If you have reached your calorie goal (by taking enough steps), a "happy face" icon will be shown. If you have not reached your calorie goal, you will see a "sad face" icon.

**Note:** While the included software can store the information of unlimited users, the pedometer can store only one user's information at a time.

**Note:** If the pedometer's buttons are not pressed for 30 seconds while adjusting your settings, it will exit Time Mode. To re-enter Time Mode, simply start from Step #1 again.

**RUNNING:** If you plan to run with the pedometer, you will need to add 12" (30.5 cm) to your Height, using the above setup instructions. The pedometer is calibrated for the stride of a walking person, so this height adjustment will ensure a more accurate step count. Once you are done running, remember to return the Height setting back to your actual height.

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

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES. OPERATION IS SUBJECT TO THE FOLLOWING TWO CONDITIONS: (1) THIS DEVICE MAY NOT CAUSE HARMFUL INTERFERENCE, AND (2) THIS DEVICE MUST ACCEPT ANY INTERFERENCE RECEIVED, INCLUDING INTERFERENCE THAT MAY CAUSE UNDESIRABLE OPERATION.

### USING THE PODEMETER

- Make sure you have activated your pedometer (by removing the clear tape from the battery compartment on the back).
- Make sure you have set your personal data. See the SET YOUR PERSONAL DATA section for more information.
- Place the pedometer somewhere on you (in your pocket, around your neck or wrist, in your hand, in your handbag, etc.).
- The pedometer will use your first ten steps of the session to calibrate itself to your stride. (These steps will not be counted.)
- Start walking/running! (**Note:** *If you are planning on running, see the RUNNING section on the previous page to set your height properly.*)
- When you are finished, you can see your data (and your data for the last seven days). See the VIEWING YOUR DATA section for more information. **To use the included software to view and track your data, insert the included Mini CD and follow the on-screen instructions to install it.**

### VIEWING YOUR DATA

#### 7-DAY RECORD CHECK

The pedometer can store 30 days of readings and display the readings of the last seven days. Follow these steps to view the readings for your steps, calories, calorie goal, distance, and stride frequency.

#### Steps (D):

- Press MODE to toggle through the available modes until Step Mode is shown. ("Step" will appear.)
- Press REC to toggle through Days 1-7. Each day's reading will be displayed with the day number.
- Press MODE to exit and return to Step Mode.

#### Calories (E):

- Press MODE to toggle through the available modes until Step Mode is shown. ("Step" will appear.)
- Press UP. "KCAL" will appear.
- Press REC to toggle through Days 1-7. Each day's reading will be displayed with the day number.
- Press UP to return to the previous "KCAL" screen. Press MODE to exit and toggle through other available modes.

**Note:** "KCAL" stands for "kilocalories," which is actually the same as "food calories" (the calories you monitor in your diet).

#### Target Kcal (F):

- Press MODE to toggle through the available modes until Step Mode is shown. ("Step" will appear.)
- Press UP twice. "TARGET" and "KCAL" will appear. The reading here is your calorie goal.
- Press MODE twice – once to exit the "KCAL" screen, once more to return to the Mode screen.

#### Distance (G):

- Press MODE to toggle through the available modes until Step Mode is shown. ("Step" will appear.)
- Press UP three times. "ML" or "KM" (whichever unit of measurement the pedometer is set to) will appear.
- Press REC to toggle through Days 1-7. Each day's reading will be displayed with the day number.
- Press MODE twice – once to exit the screen, once more to return to the Mode screen.

#### Stride Frequency (H):

- Press MODE to toggle through the available modes until Steps per Minute Mode is shown. ("S/MIN" will appear.) This number is your average steps per minute (stride frequency) across the week.
- Press MODE to exit and toggle through the other available modes.

### RESET THE PODEMETER

The pedometer's readings (steps, calories, and distance) can be reset one of two ways:

- Press MODE to toggle through the available modes until Step Mode is shown. ("Step" will appear.) Hold UP for three seconds to reset it.
- At midnight of each day (according to the pedometer's clock), the pedometer's readings will be reset automatically.

### REPLACING THE BATTERY

The included CR2032 Li battery can last four months or up to one year. To replace it:

- Use a small Phillips-head screwdriver (not included) to remove the battery compartment cover.
- Replace the old battery with a fresh CR2032 Li battery. Make sure the battery's flat side (with writing) is on the outside, facing you.
- Replace the compartment cover and screw, tightening it with your screwdriver.

**Note:** After replacing the battery, you may need to set the time again, though the pedometer should retain its other data.

### WARRANTY & SUPPORT INFORMATION

Warranty and technical support information can be found online at [www.ionhealthonline.com](http://www.ionhealthonline.com) or by calling 01252 896 080 (U.K.).

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